

Mediox

Mediox




MEDIOX Orvosi Műszergyártó Kft

3324 Felsőtárkány, 2473 Hrsz. HUNGARY
Fax: 0036 36 431 132
Tel: 0036 36 431 132

Mediox

Posterior Pedicle screw system-CTS 3.5

MEDIOX Orvosi Műszergyártó Kft




Multi-axial Screw

- Allows for 60 degree rotating, or 30 degree in any direction;
- Rotating saddle reduces rod contouring;



Low profile crosslink



Occipital plate

- Two occipital plate options available for surgeon.
- The screw is screwed into the screw hole on the body of the occipital brace and can be rotated at an angle of ±8 degrees.



Titanium Rod

- Φ3.5 titanium rod with reliable strength;
- Easily contours to meet individual anatomy;



Lateral Connector

- Easy to place of non-linear screws;
- Allow for increased angle of screw insertion;



Set Screw

- Allows for independent placement and easy surgery;
- Internal set screw optimizes placement and visualization;
- There are twisted threads and a smooth outer surface.



Parallel rod connector

Occipital connector

Hooks

Cancellous Multi Axial Screws

Product Code	D	L
950603508	Φ3.5	8
950603510	Φ3.5	10
950603512	Φ3.5	12
950603514	Φ3.5	14
950603516	Φ3.5	16
950603518	Φ3.5	18
950603520	Φ3.5	20
950603522	Φ3.5	22
950603524	Φ3.5	24
950604008	Φ4.0	8
950604010	Φ4.0	10
950604012	Φ4.0	12
950604014	Φ4.0	14
950604016	Φ4.0	16
950604018	Φ4.0	18
950604020	Φ4.0	20
950604022	Φ4.0	22
950604024	Φ4.0	24



Cortex Multi Axial Screws

Product Code	D	L
950614020	Φ4.0	20
950614022	Φ4.0	22
950614024	Φ4.0	24
950614026	Φ4.0	26
950614028	Φ4.0	28
950614030	Φ4.0	30
950614032	Φ4.0	32
950614034	Φ4.0	34
950614036	Φ4.0	36
950614038	Φ4.0	38
950614040	Φ4.0	40
950614042	Φ4.0	42
950614044	Φ4.0	44
950614046	Φ4.0	46
950614048	Φ4.0	48
950614050	Φ4.0	50
950614052	Φ4.0	52



Cortex Occipital Screws

Product Code	D	L
950633506	Φ3.5	6
950633508	Φ3.5	8
950633510	Φ3.5	10
950633512	Φ3.5	12
950633514	Φ3.5	14
950633516	Φ3.5	16
950633518	Φ3.5	18
950634006	Φ4.0	6
950634008	Φ4.0	8
950634010	Φ4.0	10
950634012	Φ4.0	12
950634014	Φ4.0	14
950634016	Φ4.0	16
950634018	Φ4.0	18
950634020	Φ4.0	20



Cortex Shaft Screws

Product Code	D	L
950623526	Φ3.5	26
950623528	Φ3.5	28
950623530	Φ3.5	30
950623532	Φ3.5	32
950623534	Φ3.5	34
950623536	Φ3.5	36
950623538	Φ3.5	38
950623540	Φ3.5	40



Cross Link

Product Code	L
950695040	40
950695045	45
950695050	50
950695055	55
950695060	60
950695065	65
950695070	70
950695075	75



Cross Link-C

Product Code	L
950695140	40
950695145	45
950695150	50
950695155	55
950695160	60
950695165	65
950695170	70
950695175	75



Rod

Product Code	D	L
950690000	Φ3.5	80
950690001	Φ3.5	120
950690002	Φ3.5	240
950690100	Φ3.5/Φ5.5	120/170
950690101	Φ3.5/Φ5.5	240/250
950690102	Φ3.5/Φ5.5	190/500

Occipital Plates

Product Code	Specification
950691000	3 holes
950691100	4 holes



Hooks

Product Code	Specification
950692000	Small
950692100	Large



Set Screw

Product Code	D
950694000	Φ6



Occipital Connector

Product Code	Specification
950693300	



Open Transverse Connector

Product Code	Specification
950693000	



Close Transverse Connector

Product Code	Specification
950693100	



Parallel Rod Connector

Product Code	Specification
950693200	

